

Product Datasheet

Mouse Src Antibody (N-term) (orb1934137)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Src antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 98-125 amino acids from the N-terminal region of mouse Src.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59891 Da
Uniprot ID	P05480
NCBI	NP_033297.2, NP_001020566.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.