

Product Datasheet

Mouse Hoxc4 Antibody (Center) (orb1934141)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Hoxc4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 139-166 amino acids from the Central region of mouse Hoxc4.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	29865 Da
Uniprot ID	Q08624
NCBI	NP_038581.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.