

Product Datasheet

AP4B1 Antibody (Center) (orb1934166)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This AP4B1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 492-520 amino acids from the Central region of human AP4B1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83260 Da
Uniprot ID	Q9Y6B7
NCBI	NP_006585.2 , NP_001240782.1 , NP_001240781.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.