

Product Datasheet

Fbn2 antibody (orb19345)

Description

Goat polyclonal antibody to Fbn2

Species/Host

Goat

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

Fibrillin 2 (mouse)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 314kDa according to NP_034311.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

KELRKLEEHNEDD

Clonality

Polyclonal

MW

314

Dilution Range

ELISA: 1:2000

Expiration Date

12 months from date of receipt.