

Product Datasheet

KLRG2 Antibody (N-term) (orb1934984)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This KLRG2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 20-49 amino acids from the N-terminal region of human KLRG2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42851 Da
Uniprot ID	A4D1S0
NCBI	NP_940910.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.