

Product Datasheet

Mouse Adrbk2 Antibody (C-term) (orb1935462)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This Mouse Adrbk2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 524-552 amino acids from the C-terminal region of mouse Adrbk2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	79657 Da
Uniprot ID	Q3UYH7
NCBI	NP_796052.2
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.