

Product Datasheet

SLC22A16 antibody (orb19355)

Description

Goat polyclonal antibody to SLC22A16

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

SLC22A16

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments gave an approx 85-90kDa band in Human Testis lysates after 0.5 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 64.6kDa according to NP_149116.2. The 85-90kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

CSRNKRENTSSL

Clonality

Polyclonal

MW


64.6

Entrez
85413
Dilution Range

ELISA: 1:128000, IHC-P: 1-2µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Testis