

Product Datasheet

RCAN3 Antibody (C-term) (orb1935844)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This RCAN3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 194-222 amino acids from the C-terminal region of human RCAN3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27492 Da
Uniprot ID	Q9UKA8
NCBI	NP_001238911.1, NP_001238912.1, NP_038469.1, NP_001238908.1, NP_001238909.1, NP_001238907.1, NP_001238910.1, NP_001238914.1, NP_001238906.1
Dilution Range	WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.