

Product Datasheet

Mouse Jak3 Antibody (N-term) (orb1935976)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This Mouse Jak3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 133-161 amino acids from the N-terminal region of mouse Jak3.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	122641 Da
Uniprot ID	Q62137
NCBI	NP_001177759.1 , NP_034719.2
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.