


Product Datasheet

GATA5 antibody (orb19361)

Description	Goat polyclonal antibody to GATA5
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GATA5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band observed in lysates of liver cell line HepG2 (calculated MW of 41.3kDa according to NP_536721.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	ELAARPGWAQTATAD
Clonality	Polyclonal
MW	41.3
Entrez	140628
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...