

Product Datasheet

NFIA Antibody (C-term) (orb1936104)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	This NFIA antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 470-498 amino acids from the C-terminal region of human NFIA.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55944 Da
Uniprot ID	Q12857
NCBI	NP_005586.1 , NP_001138983.1 , NP_001128145.1 , NP_001138984.1
Dilution Range	IF: 1:10~50, WB: 1:1000
Expiration Date	12 months from date of receipt.