

Product Datasheet

TRPV1 Antibody (N-term) (orb1936238)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This TRPV1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 124-153 amino acids from the N-terminal region of human TRPV1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	94956 Da
Uniprot ID	Q8NER1
NCBI	NP_542437.2, NP_061197.4, NP_542435.2, NP_542436.2
Dilution Range	WB: 1:2000, WB: 1:1000, IHC-P-Leica: 1:500, IHC-P-Leica: 1:500
Expiration Date	12 months from date of receipt.