

Product Datasheet

ADNP2 Antibody (Center) (orb1936243)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This ADNP2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 702-731 amino acids from the Central region of human ADNP2.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	122833 Da
Uniprot ID	Q6IQ32
NCBI	NP_055728.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.