

Product Datasheet

Bad BH3 Domain Antibody (orb1936829)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	This Bad BH3 Domain antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 92-127 amino acids from human Bad BH3 Domain.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18392 Da
Uniprot ID	Q92934
NCBI	NP_116784.1 , NP_004313.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.