

Product Datasheet

CACNB4 antibody (orb19370)

Description	Goat polyclonal antibody to CACNB4
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CACNB4 (internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 55-60kDa band observed in Human Bone Marrow lysates (calculated MW of 54.7kDa according to NP_001005747.1 and 58.2kDa according to NP_000717.2). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	DYPDSYQDTYKPH
Clonality	Polyclonal
MW	54.7; 58.2; 56.4
Entrez	785
Dilution Range	ELISA: 1:4000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



orb19370
(0.2ug/ml)
staining of
Human Bo...