

Product Datasheet

TAC4 Antibody (C-term) (orb1937143)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This TAC4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 85-113 amino acids from the C-terminal region of human TAC4.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	12305 Da
Uniprot ID	Q86UU9
NCBI	NP_733786.2, NP_001070971.1, NP_001070973.1, NP_001070972.1, NP_001070974.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.