

Product Datasheet

K0746 Antibody (N-term) (orb1937170)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	This K0746 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 253-282 amino acids from the N-terminal region of human K0746.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	128567 Da
Uniprot ID	Q68CR1
NCBI	NP_056002.2
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.