

Product Datasheet

CDRT15L1 Antibody (Center) (orb1937387)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This CDRT15L1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 134-162 amino acids from the Central region of human CDRT15L1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30100 Da
Uniprot ID	A8MXV6
NCBI	NP_001177719.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.