

Product Datasheet

PSAPL1 Antibody (N-term) (orb1937503)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	This PSAPL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 130-158 amino acids from the N-terminal region of human PSAPL1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56627 Da
Uniprot ID	Q6NUJ1
NCBI	NP_001078851.1
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.