

Product Datasheet

C109B Antibody (C-term) (orb1937540)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This C109B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 291-320 amino acids from the C-terminal region of human C109B.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39082 Da
Uniprot ID	Q9NWR8
NCBI	NP_060388.2
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.