

## Product Datasheet

### PCDHB3 Antibody (N-term) (orb1937702)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Affinity Purified Rabbit Polyclonal Antibody (Pab)                                                                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IF, IHC-P, WB                                                                                                                                                                |
| <b>Target</b>              | This PCDHB3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 12-39 amino acids from the N-terminal region of human PCDHB3.       |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 86714 Da                                                                                                                                                                     |
| <b>Uniprot ID</b>          | <b>Q9Y5E6</b>                                                                                                                                                                |
| <b>NCBI</b>                | <b>NP_061760.1</b>                                                                                                                                                           |
| <b>Dilution Range</b>      | IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |