

Product Datasheet

COL21A1 Antibody (N-term) (orb1937708)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This COL21A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 204-232 amino acids from the N-terminal region of human COL21A1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	99369 Da
Uniprot ID	Q96P44
NCBI	NP_110447.2
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.