

Product Datasheet

PCDHB13 Antibody (Center) (orb1937790)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This PCDHB13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 273-300 amino acids from the Central region of human PCDHB13.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87552 Da
Uniprot ID	Q9Y5F0
NCBI	NP_061756.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.