

Product Datasheet

PCDHGA8 Antibody (C-term) (orb1937794)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This PCDHGA8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 781-808 amino acids from the C-terminal region of human PCDHGA8.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	101480 Da
Uniprot ID	Q9Y5G5
NCBI	NP_114477.1 , NP_054723.1
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.