

Product Datasheet

PPP1R3G Antibody (C-term) (orb1937810)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	This PPP1R3G antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 269-298 amino acids from the C-terminal region of human PPP1R3G.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38019 Da
Uniprot ID	B7ZBB8
NCBI	NP_001138587.1
Dilution Range	WB: 1:1000, FC: 1:10~50
Expiration Date	12 months from date of receipt.