

Product Datasheet

WDR42C Antibody (N-term) (orb1937879)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This WDR42C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 112-140 amino acids from the N-terminal region of human WDR42C.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	71191 Da
Uniprot ID	POC7V8
NCBI	NP_001130005.1
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.