

Product Datasheet

ARL5C Antibody (N-term) (orb1937940)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This ARL5C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 34-61 amino acids from the N-terminal region of human ARL5C.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20591 Da
Uniprot ID	A6NH57
NCBI	NP_001137440.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.