

Product Datasheet

GOT1L1 Antibody (C-term) (orb1937943)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This GOT1L1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 346-375 amino acids from the C-terminal region of human GOT1L1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47305 Da
Uniprot ID	Q8NHS2
NCBI	NP_689626.2
Dilution Range	IF: 1:10~50, WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.