

Product Datasheet

C6orf138 Antibody (Center) (orb1938222)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This C6orf138 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 563-591 amino acids from the Central region of human C6orf138.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	96371 Da
Uniprot ID	Q6ZW05
NCBI	NP_001013754.3, NP_997382.2
Dilution Range	WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.