

Product Datasheet

LOC441476 Antibody (Center) (orb1938224)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This LOC441476 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 134-162 amino acids from the Central region of human LOC441476.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42496 Da
Uniprot ID	Q8N7X2
NCBI	NP_001243628.1, NP_001004353.2, NP_001243630.1, NP_001243629.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.