


Product Datasheet

GNRHR antibody (orb19383)

Description	Goat polyclonal antibody to GNRHR
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GRHR / LRHR
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in Human lysates of Peripheral Blood Mononucleocytes (calculated MW of 37.7kDa according to NP_000397.1). The observed molecular weight corresponds to earlier findings with products from other commercial sources. Recommended concentration: 0.5-1.5 µg/ml.
Protein Sequence	QKWTQKKEKGK
Clonality	Polyclonal
MW	37.7; 27.7
Entrez	2798
Dilution Range	ELISA: 1:32000, WB: 0.5-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb19383
(0.5ug/ml)
staining of
Human Pe...