

Product Datasheet

FAM108C1 Antibody (Center) (orb1938302)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This FAM108C1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 92-121 amino acids from the Central region of human FAM108C1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35831 Da
Uniprot ID	Q6PCB6
NCBI	NP_067037.1
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.