

Product Datasheet

LENG9 Antibody (N-term) (orb1938321)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	This LENG9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 127-155 amino acids from the N-terminal region of human LENG9.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53167 Da
Uniprot ID	Q96B70
NCBI	NP_945339.2
Dilution Range	WB: 1:1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.