

Product Datasheet

CCDC126 Antibody (C-term) (orb1938652)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	This CCDC126 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 77-106 amino acids from the C-terminal region of human CCDC126.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15668 Da
Uniprot ID	Q96EE4
NCBI	NP_620126.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.