

Product Datasheet

WWC1 Antibody (N-term) (orb1938741)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This WWC1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 286-313 amino acids from the N-terminal region of human WWC1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	125301 Da
Uniprot ID	Q8IX03
NCBI	NP_001155134.1, NP_056053.1, NP_001155133.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.