

Product Datasheet

OSBPL3 antibody (orb19389)

Description	Goat polyclonal antibody to OSBPL3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	ORP3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.
Protein Sequence	MSDEKNLGVSQKL
Clonality	Polyclonal
MW	101; 97.7; 97.3; 93.8
Entrez	26031
Dilution Range	ELISA: 1:32000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.