



Product Datasheet

NBEA antibody (orb19390)

Description	Goat polyclonal antibody to NBEA
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Neurobeachin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 300kDa band observed in lysates of transfected HEK293 and endogenously in lysates of insulinoma cell line B-TC3 (calculated MW of 327kDa according to NP_056493.3). Data kindly provided by Professor John Creemers, K.U. Leuven, Belgium. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DFNRWHYEHQNRY
Clonality	Polyclonal
MW	328; 84.0
Entrez	26960
Dilution Range	ELISA: 1:4000, IHC-P: 5-10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells u...



Immunohistochemical staining of Human Ce...