

Product Datasheet

HNF1A antibody (orb19391)

Description	Goat polyclonal antibody to HNF1A
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	TCF1 / HNF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Preliminary experiments in human liver and kidney lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	DELPTKKGRRNRFK
Clonality	Polyclonal
MW	67.4
Entrez	6927
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml).
Expiration Date	12 months from date of receipt.