

Product Datasheet

SERPINC1 Antibody (orb1939271)

Description	Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	This SERPINC1 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 364-393 amino acids from human SERPINC1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52602 Da
Uniprot ID	P01008
NCBI	NP_000479.1
Dilution Range	WB: 1:500~1000
Expiration Date	12 months from date of receipt.