

Product Datasheet

NUCL Antibody (orb1939424)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	This NUCL Monoclonal antibody is generated from mice immunized with a protein from human NUCL.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1?
Clonality	Monoclonal
MW	76614 Da
Uniprot ID	P19338
NCBI	NP_005372.2
Dilution Range	IF: 1:10~50, WB: 1:500~1000, IHC-P: 1:50~100
Expiration Date	12 months from date of receipt.