

Product Datasheet

CT45A4/CT45A3/CT45A2/CT45A1/CT45A6 Antibody (Center) (orb1939634)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This CT45A4/CT45A3/CT45A2/CT45A1/CT45A6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 106-135 amino acids from the Central region of human CT45A4/CT45A3/CT45A2/CT45A1/CT45A6.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q8N7B7
NCBI	NP_001017436.1, NP_001017438.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.