

Product Datasheet

Beta Synuclein Pre-formed Fibrils (Type 1) (orb1939724)

Description	Human Recombinant Beta Synuclein Protein Pre-formed Fibrils (Type 1)
Reactivity	Human
Tested Applications	In vitro, In vivo, SDS-PAGE, WB
Target	Beta Synuclein
Preservatives	PBS pH 7.4
Concentration	Lot/batch specific. See included datasheet.
Storage	-80°C
Tag	No tag
Note	For research use only
Application notes	Certified > 95% pure using SDS-PAGE analysis.
Purity	> 95%
Source	Recombinant
MW	14.28 kDa
Uniprot ID	Q16143
Expiration Date	6 months from date of receipt.