

Product Datasheet

Anti-RILPL1 Antibody (orb1939860)

Description
Biorbyt Anti-RILPL1 Antibody catalog #

	This antibody react...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E.coli-derived human RILPL1 recombinant protein (Position: E3-K324). Human RILPL1 shares 93.2% and 93.5% amino acid (aa) sequence identity with mouse and rat RILPL1, respectively.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	Q5EBL4
Dilution Range	Western blot, 0.25-0.5 µg/ml, Human, Rat ELISA, 0.1-0.5 µg/ml, Human
Expiration Date	12 months from date of receipt.