

Product Datasheet

Anti-Bim/Bcl2l11 Antibody (orb1939869)

Description

Biorbyt Anti-Bim/Bcl2l11 Antibody catalog

	FCM(Intracellular), ELISA applic...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	E.coli-derived mouse Bim/Bcl2l11 recombinant protein (Position: Q18-R186). Mouse Bcl2l11 shares 86.7% and 98.8% amino acid (aa) sequence identity with human and rat Bcl2l11, respectively.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	054918
Dilution Range	Western blot, 0.25-0.5 µg/ml, Mouse Flow Cytometry (Fixed), 1-3 µg/1x10 ⁶ cells, Mouse ELISA, 0.1-0.5 µg/ml, Mouse
Expiration Date	12 months from date of receipt.