


Product Datasheet

AGTR1 antibody (orb19401)

Description	Goat polyclonal antibody to AGTR1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AGTR1 / AT1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 38kDa band observed in Liver lysates of Human, Mouse and Rat (calculated MW of 41.1kDa according to NP_000676.1, NP_004826.2, NP_033611.1, NP_114038.1 and NP_114438.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	EIQKNKPRNDDIFK
Clonality	Polyclonal
MW	41.1
Entrez	185
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Liver lys...