

## Product Datasheet

### CRHR1 antibody (orb19405)

## Description

Goat polyclonal antibody to CRHR1

## Species/Host

Goat

## Reactivity

Bovine, Canine, Human, Mouse, Rat, Sheep

## Conjugation

Unconjugated

## Tested

ELISA, IF, IHC, WB

## Applications

## Target

CRHR1 / CRF-R (aa 250 - 263)

## Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

## Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

## Note

For research use only

## Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: This product is successfully used in immunohistochemistry on mouse brain.WB: Preliminary experiments in human colon, duodenum and placenta lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysatesExperiment Notes:  
Immunofluorescence: This product is successfully used in immunofluorescence on mouse brain.

## Protein Sequence

KLYYDNEK CWFGKR

## Clonality

Polyclonal

## MW

47.7

## Entrez


**1394**

## Dilution Range


ELISA: 1:1000, IHC-P: 5 µg/ml

## Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Co...



Immunofluorescence analysis of MCF7 cell...