

Product Datasheet

DISC1 antibody (orb19406)

Description	Goat polyclonal antibody to DISC1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DISC1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 100kDa band observed in human brain (cerebellum, hippocampus and frontal cortex) lysates (calculated MW of 93.6kDa according to NP_001012975.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HSEDKRKTPCLKVLE
Clonality	Polyclonal
MW	93.6; 91.1
Entrez	27185
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19406
(0.1ug/ml)
staining of
human br...