


Product Datasheet

TBC1D4 antibody (orb19407)

Description	Goat polyclonal antibody to TBC1D4
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AS160 / TBC1D4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Skeletal Muscle shows striated pattern staining of majority of muscular fibres. Recommended concentration, 3-5ug/ml.WB: Approx 150kDa band observed in lysates of cell lines Daudi and MOLT4 (calculated MW of 147kDa according to NP_055647.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DDPEKIEERKKSK
Clonality	Polyclonal
MW	147
Entrez	9882
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



orb19407
(0.1ug/ml)
staining of
Daudi ce...



orb19407
(3.8ug/ml)
staining of
paraffin...