

Product Datasheet

FADS1 antibody (orb19411)

Description	Goat polyclonal antibody to FADS1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FADS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 50kDa band observed in Human Liver and Lung lysates (calculated MW of 51.9kDa according to NP_037534.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QPSFEPTKNKELTDE
Clonality	Polyclonal
MW	57.8
Entrez	3992
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung lysa...



Western blot analysis of Human Lung lysa...