


Product Datasheet

TCF3 antibody (orb19418)

Description	Goat polyclonal antibody to TCF3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TCF3 / ITF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Kidney shows strong nuclear staining in DCT cells. Recommended concentration, 3-5ug/ml.WB: Approx 75kDa band observed in lysates of cell line Daudi (calculated MW of 67.6kDa according to NP_003191.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	KAPRARTSPDEDED
Clonality	Polyclonal
MW	67.6
Entrez	6929
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



orb19418
(0.01ug/ml)
staining of
Daudi c...



orb19418
(3.8ug/ml)
staining of
paraffin...